


# HPLC REPORT

Sample: K-15 Analyzed date: 2014-03-31  
 Analyst: HXH  
 Lot. No.: P140302-WF566269  
 Column: SinoChrom ODS-BP, 4.6\*250mm, 5µm  
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
 Gradient:
 

	A	B
0.0min	11%	89%
25.0min	36%	64%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl  
 Wavelength: 220nm  
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	10.067	143617.234	3.4635
2	10.590	3661968.000	88.3116
3	10.590	188157.719	4.5376
4	10.998	145011.656	3.4971
5	16.207	7890.756	0.1903
Total		4146645.365	100